DOC022.98.80204

2100N IS

05/2015, Edition 6

Basic User Manual
Grundlegende Bedienungsanleitung
Manuale di base per l'utente
Manuel d'utilisation de base
Manual básico del usuario
Manual de operações básicas
Základní návod k použití
Grundlæggende brugerhåndbog
Basisgebruikershandleiding
Podstawowa instrukcja obsługi
Grundläggande bruksanvisning
Peruskäyttöohje
Manualul de bază al utilizatorului
Basic User Manual
Βασικό εγχειρίδιο λειτουργίας

English.....	3
Deutsch.....	17
Italiano.....	32
Français.....	47
Español.....	62
Português.....	77
Čeština.....	92
Dansk.....	106
Nederlands.....	120
Polski.....	135
Svenska.....	149
Suomi.....	163
Română.....	177
Türkçe.....	191
Ελληνικά.....	205

Table of contents

[Specifications](#) on page 3 [Standard operation](#) on page 9
[General information](#) on page 4 [Maintenance](#) on page 14
[User interface](#) on page 7 [Troubleshooting](#) on page 14
[Startup](#) on page 8

Specifications

Specifications are subject to change without notice.

Specification	Details
Measurement method	Nephelometric
Regulatory	Meets ISO 7027, DIN EN 27027, DIN 38404 and NFT 9033 ASTM D7315 - Standard Test Method for Determination of Turbidity Above 1 Turbidity Unit (TU) in Static Mode ASTM D6655 - Standard Test Method for Determination of Turbidity Below 5 NTU in Static Mode
Light source	Light-emitting diode (LED) at 860 ± 30 nm
Measurement modes	FNU and NTU
Range	FNU (manual range): 0–0.999, 0–9.99, 0–99.9, 0–1000 FNU (auto range): 0–1000 NTU (manual range): 0–0.999, 0–9.99, 0–99.9, 0–1000 NTU (auto range): 0–1000
Accuracy ^{1, 2}	±2% of reading plus 0.01 FNU/NTU from 0–1000 FNU/NTU

¹ Turbidity specifications identified using recently prepared formazin standard and matched 1-inch sample cells.

² Intermittent electromagnetic radiation of 3 volts/meter or greater may cause slight accuracy shifts.

Specification	Details
Resolution	Turbidity: 0.001 FNU/NTU (on lowest range)
Repeatability	±1% of reading or 0.01 FNU, whichever is greater (under reference conditions)
Response time	Signal averaging off: 6.8 seconds Signal averaging on: 14 seconds (when 10 measurements are used to calculate the average)
Stabilization time	Immediately
Reading modes	Manual or auto range, signal averaging on or off
Power requirement	115–230 VAC, 50/60 Hz (automatic power selection) 28 W maximum
Pollution degree/installation category	2; II
Protection Class	1
Operating conditions	Temperature: 0 to 40 °C (32 to 104 °F) Relative humidity: 0–90% at 25 °C, 0–75% at 40 °C, noncondensing Altitude: 2000 m (6560 ft) maximum Indoor use only
Storage conditions	–40 to 60 °C (–40 to 140 °F), instrument only
Interface	RS232C serial interface by way of DB9 subminiature D-shell connector for data output to computer or printer, and data input (command). No handshaking.

Specification	Details
Air purge	Dry nitrogen or instrument grade air (ANSI MC 11.1, 1975) 0.1 scfm at 69 kPa (10 psig); 138 kPa (20 psig) maximum Hose barb connection for 1/8-inch tubing
Sample cells	Round cells 95 x 25 mm (3.74 x 1 in.) borosilicate glass with rubber-lined screw caps Note: Smaller sample cells (less than 25 mm) can be used when a cell adapter is used.
Sample requirements	25 mm sample cell: 20 mL minimum 0 to 95 °C (32 to 203 °F)
Enclosure	High-impact polycarbonate plastic
Dimensions	30.5 x 40 x 15.6 cm (12.0 x 15.7 x 6.1 in.)
Weight	3.8 kg (8.5 lb)
Certification	CE, cETLus

General information

In no event will the manufacturer be liable for direct, indirect, special, incidental or consequential damages resulting from any defect or omission in this manual. The manufacturer reserves the right to make changes in this manual and the products it describes at any time, without notice or obligation. Revised editions are found on the manufacturer's website.

Additional information

Additional information is available on the manufacturer's website.

Safety information

NOTICE

The manufacturer is not responsible for any damages due to misapplication or misuse of this product including, without limitation, direct, incidental and consequential damages, and disclaims such damages to the full extent permitted under applicable law. The user is solely responsible to identify critical application risks and install appropriate mechanisms to protect processes during a possible equipment malfunction.

Please read this entire manual before unpacking, setting up or operating this equipment. Pay attention to all danger and caution statements. Failure to do so could result in serious injury to the operator or damage to the equipment.

Make sure that the protection provided by this equipment is not impaired. Do not use or install this equipment in any manner other than that specified in this manual.

Use of hazard information

⚠ DANGER

Indicates a potentially or imminently hazardous situation which, if not avoided, will result in death or serious injury.

⚠ WARNING

Indicates a potentially or imminently hazardous situation which, if not avoided, could result in death or serious injury.

⚠ CAUTION

Indicates a potentially hazardous situation that may result in minor or moderate injury.

NOTICE

Indicates a situation which, if not avoided, may cause damage to the instrument. Information that requires special emphasis.

Precautionary labels

Read all labels and tags attached to the instrument. Personal injury or damage to the instrument could occur if not observed. A symbol, if noted

on the instrument, will be included with a danger or caution statement in the manual.

	This symbol, if noted on the instrument, references the instruction manual for operation and/or safety information.
	Electrical equipment marked with this symbol may not be disposed of in European public disposal systems after 12 August of 2005. In conformity with European local and national regulations (EU Directive 2002/96/EC), European electrical equipment users must now return old or end-of-life equipment to the Producer for disposal at no charge to the user. Note: For return for recycling, please contact the equipment producer or supplier for instructions on how to return end-of-life equipment, producer-supplied electrical accessories, and all auxiliary items for proper disposal.

Certification

Canadian Radio Interference-Causing Equipment Regulation, IECIS-003, Class A:

Supporting test records reside with the manufacturer.

This Class A digital apparatus meets all requirements of the Canadian Interference-Causing Equipment Regulations.

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FCC Part 15, Class "A" Limits

Supporting test records reside with the manufacturer. The device complies with Part 15 of the FCC Rules. Operation is subject to the following conditions:

1. The equipment may not cause harmful interference.
2. The equipment must accept any interference received, including interference that may cause undesired operation.

Changes or modifications to this equipment not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment. This equipment has been tested and found to comply with the limits for a Class A digital device, pursuant to Part 15 of the FCC rules. These limits are designed to provide reasonable

protection against harmful interference when the equipment is operated in a commercial environment. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instruction manual, may cause harmful interference to radio communications. Operation of this equipment in a residential area is likely to cause harmful interference, in which case the user will be required to correct the interference at their expense. The following techniques can be used to reduce interference problems:


1. Disconnect the equipment from its power source to verify that it is or is not the source of the interference.
2. If the equipment is connected to the same outlet as the device experiencing interference, connect the equipment to a different outlet.
3. Move the equipment away from the device receiving the interference.
4. Reposition the receiving antenna for the device receiving the interference.
5. Try combinations of the above.

Product overview

The 2100N IS laboratory turbidimeter measures turbidity in FNU (Formazin nephelometric units) and NTU (nephelometric turbidity units). NTUs are calculated using the conversion factors of 1.0 NTU per 1.0 FNU.


The turbidimeter has an RS232 output for connection to a printer, data logger or computer.

Figure 1 Front overview

1 Keypad	4 Cover for the sample cell compartment
2 Sample cell holder	5 Five-digit LED display
3 Light shield	

Figure 2 Back overview


1 Power cord connector	4 DB9 connector for RS232 cable
2 Fuse holder	5 Air purge fitting
3 Power switch	

Product components

Refer to [Figure 3](#) to make sure that all components have been received. If any of these items are missing or damaged, contact the manufacturer or a sales representative immediately.

Figure 3 Instrument components


1 2100N IS turbidimeter	5 StablCal® Calibration kit
2 Oiling cloth	6 Gelex® secondary turbidity standardization kit ³
3 Six 1" sample cells (30 mL) with caps	7 Dust cover
4 Silicone oil	8 Power cord

³ Supplied with 4790000 only.

User interface

Figure 4 Keypad

1 ENTER key	5 UNITS AVG key key
2 EDIT (arrow) keys	6 SETUP key
3 RANGE key	7 PRINT key
4 UNITS/Exit key	8 CAL key


Table 1 Key descriptions

Key	Description
ENTER	Enters the value on the display. Starts the measurement of a calibration standard. Clears data from the buffer.
EDIT (arrow keys)	Changes the numbers and/or letters on the display. Steps through the calibration standards. The right arrow key moves the cursor to the previous or next digit.
RANGE	Selects automatic or manual ranging.
UNITS Exit	Selects the unit of measure. Exits Calibration or Setup mode without saving changes.

Table 1 Key descriptions (continued)

Key	Description
UNITS AVG	Turns signal averaging on or off.
SETUP	Turns on Setup mode and starts the selection of the setup number.
PRINT	Sends the data that is on the display to a printer or computer. Sends a calibration data report to a printer or computer when in Calibration mode. Sends diagnostic results to a printer or computer if held down when the instrument is turned on. Provides a print of the setup commands when in Setup mode.
CAL	Starts or completes a calibration.

Figure 5 Indicator lights

1 NTU light	6 CAL? light
2 FNU light	7 UNITS AVG light
3 Lamp light	8 SETUP light
4 Manual Range light	9 S0–S3 lights
5 Auto Range light	

Table 2 Light descriptions

Light	Description
NTU	Illuminated when the instrument is set for NTU unit of measure.
FNU	Illuminated when the instrument is set for FNU unit of measure.
	Illuminated when the instrument light source is on. Flashes when there is not sufficient light for measurement.

Table 2 Light descriptions (continued)

Light	Description
Manual Range	Illuminated when the instrument is in manual ranging mode.
Auto Range	Illuminated when the instrument is in auto ranging mode.
CAL?	Turns on during a calibration if the calibration data is not within the acceptable range. Flashes when the instrument should be calibrated. Note: The CAL? light applies when a 25-mm sample cell is used. Ignore the CAL? light if illuminated during calibration when a smaller sample cell is used. Push UNITS/Exit to start measurements.
UNITS AVG	Illuminated when signal averaging is on.
SETUP	Illuminated when Setup mode is selected.
S0–S3	Show the current calibration point standard that is in use during calibration.

Startup

Turn the instrument on

CAUTION

Infrared Light Hazard. The infrared light produced by this instrument can cause eye injury. The infrared light source in this instrument only receives power when the sample cell cover is closed.

- Put the instrument on a stable, level surface that is free of vibration. Do not put in direct sunlight.
- Make sure that there is air circulation around the instrument. Keep the back and area below the instrument free of material that could decrease air flow through the vents.
- Connect the power cord to the power plug on the back of the instrument.

4. Connect the power cord to a power socket with ground contact.
5. Push the power switch on the back of the instrument to turn the instrument on.

Turn the keypad sound off (optional)

By default, the instrument makes an audible sound when a key is pushed. To turn the keypad sound off:

1. Push **SETUP**. The SETUP light turns on.
2. Use the arrow keys to select 00.
3. Push **ENTER**.
4. Use the arrow keys to select the sound option:

Option	Description
BP	An audible sound is made when a key is pushed.
BP OF	No sound is made when a key is pushed.

5. Push **ENTER**.
6. Push **SETUP**.

Standard operation

Calibrate the turbidimeter with StablCal® Standards

Calibrate the turbidimeter before it is used for the first time using the StablCal® sealed vial standards provided. As an alternative, calibration can be done with recently prepared formazin standards.

Calibrate the turbidimeter at least every 3 months or as specified by the regulating authority when data is used for ISO 7027 reporting.

Note: Unknown results may occur if standards other than the recommended calibration points are used. The recommended calibration points (< 0.1, 20, 200 and 1000 NTU) provide the best calibration accuracy. Use of standards other than StablCal, or user-prepared formazin, may result in less accurate calibrations. The manufacturer cannot guarantee the performance of the instrument if calibrated with co-polymer styrenedivinylbenzene beads or other suspensions.

Prepare the StablCal standards

When received and at intervals:

1. Clean the exterior surface of the StablCal vials with laboratory glass cleaning detergent.
2. Rinse the vials with distilled or deionized water.
3. Dry the vials with a lint-free cloth.

Note: Never shake or invert the < 0.1 NTU standard. If the standard has been mixed or shaken, do not move the vial for 15 minutes or more before using.

Note: Do not remove the caps from the sealed vials.

Make sure that the StablCal standards are at ambient instrument temperature before use (and no greater than 40 °C (104 °F)).

Mix the standards before use:

1. Open the case lid. Remove the < 0.1 NTU standard from the plastic case.
2. Leave the other standards in the case. Close the case lid.
3. Shake the case vigorously for at least 10 seconds.
4. Let the standards stand with no movement for 3–5 minutes before use.

Calibration notes

- Make sure that the instrument is in the same ambient conditions as where it is used.
- Make sure that the standards are at the same ambient temperature as the instrument before use.
- Use only the provided silicone oil. This silicone oil has the same refractive index as the vial glass and masks minor glass differences and scratches.
- Store the oiling cloth in a plastic storage bag to keep the cloth clean.
- If power is lost during calibration, the new calibration data is lost and the last calibration data is used. To exit a calibration and not save the new values, push **UNITS/Exit**.
- In Calibration mode, automatic range and signal averaging on are selected. When calibration is completed, all operational modes go back to the last settings.

- All nephelometric (turbidity units of measure) calibrations are done at the same time.
- The 4000-NTU standard does not have to be measured during calibration if FNU values will be measured. Push **CAL** after the 1000 NTU standard is measured to complete the calibration procedure.

- The FNU values of StablCal standards and formazin standards are calculated using the conversion factors of 1 FNU = 1 NTU.

StablCal calibration procedure


1. Push **CAL**.

The S0 light turns on. The NTU value of the dilution water used in the last calibration is shown on the display.


2. Get the < 0.1 NTU vial. Clean the vial with a soft, lint-free cloth to remove water spots and fingerprints. Do not invert the vial.

3. Apply a small bead of silicone oil from the top to the bottom of the vial.

4. Use the oiling cloth to apply the oil equally to the surface of the vial. Remove the excess oil. Make sure that the vial is almost dry.


5. Put the vial in the sample cell holder with the triangle on the vial aligned with the reference mark on the sample cell holder. Close the cover.


6. Push **ENTER**.

The instrument display counts down, then measures the standard. The next expected standard (e.g., 20.00) is shown. The S1 light turns on.

7. Remove the vial from the sample cell holder.

8. Do steps 5–10 for the other StablCal vials (from lowest to highest NTU standard).

The S0 light turns on after the last vial is measured.

9. Push **CAL**.

The instrument saves the new calibration data and goes back to Measurement mode.

StablCal standards storage

- Do not move a StablCal standard to a different container for storage. Keep StablCal standards in the plastic case provided with the cover closed.
- Store at 5 to 25 °C (41 to 77 °F).
- For long-term storage (more than one month between use), keep at 5 °C (41 °F).

Turbidity measurement

⚠ WARNING

Potential explosion and fire hazard. This instrument is for measuring water based samples. Do not measure solvent or combustible based samples.

For accurate turbidity readings use clean sample cells and remove air bubbles.

Measurement notes

Proper measurement techniques are important in minimizing the effects of instrument variation, stray light and air bubbles. For accurate and repeatable measurements:

Instrument

- Make sure that the instrument is on a level, stationary surface that is free of vibration during the measurement.
- Instrument stabilization is immediate. No warm-up time is necessary.
- Always close the sample compartment lid during measurement, calibration and storage.
- Remove the sample cell from the instrument and turn off the instrument if the instrument is stored for an extended time period (more than a month).
- Keep the sample compartment lid closed to keep dust and dirt out.

Sample cells

- Always cap the sample cell to prevent spillage of the sample into the instrument.
- Always use clean sample cells in good condition. Dirty, scratched or damaged cells can result in readings that are not accurate.
- Make sure that cold samples do not “fog” the sample cell.
- Store sample cells filled with distilled or deionized water and cap tightly.

- For the best accuracy, use a single sample cell for every measurement or a flow cell.

Note: As an alternative, matched sample cells may be used for measurements but do not provide as good of accuracy or precision as a single indexed sample cell or flow cell. When using matched sample cells, align the orientation mark on the sample cell with the reference mark on the sample cell holder.

Measurement


Turbidity measurement procedure

1. Rinse a clean, empty sample cell two times with the solution to be measured and drain to waste. Fill to the line (about 30 mL) with sample and immediately put the cap on the sample cell.


2. Clean the sample cells with a soft, lint-free cloth to remove water spots and fingerprints.


3. Apply a small bead of silicone oil from the top to the bottom of the sample cells.


4. Use the oiling cloth provided to apply the oil equally to the surface of the sample cells. Remove the excess oil. Make sure that the sample cells are almost dry.

5. Gently and slowly invert the sample cell to fully mix the sample. Be careful not to add air bubbles.

6. Put the sample cell in the sample cell holder with the triangle on the sample cell aligned with the reference mark on the sample cell holder. Close the cover.

7. Read and record the value when stable.

Note: To send (via RS232) a measurement record, push **PRINT**.

Measurement techniques

Measurements may be made with different operation mode settings and optional accessories.

Calibrate the instrument whenever the sample cell pathlength is changed.

Manual or automatic ranging

The manufacturer recommends that ranging be set to automatic for most measurements.

The setting can be changed at any time during sample measurement.

Push **RANGE** repeatedly to step the instrument from automatic ranging to manual ranging and then scroll through the manual range settings.

The Manual Range light turns on when manual ranging is selected. The Auto Range light turns on when automatic ranging is selected.

Notes:

- When manual ranging is selected, the display flashes all 9s when the sample being measured is greater than the selected range. The display flashes all 0s when the sample measured is less than the selected range.
- When automatic ranging is selected, the display flashes 9s when the sample is greater than the maximum range of the instrument.

- When automatic ranging is selected, the display flashes all 0s when the measurement is less than the range of the instrument or a negative value. Calibrate the instrument.

Signal averaging on or off

Signal averaging corrects for reading fluctuations that are caused by random drifting particles in the sample. When signal averaging is on, an average reading is calculated every 3 seconds and shown on the display. The last ten measurements are used to calculate the average reading.

The manufacturer recommends that signal averaging be on for most measurements.

Push **UNITS AVG** to turn signal averaging on or off. The UNITS AVG light turns on when signal averaging is on.

Push **ENTER** when signal averaging is on to erase data in the signal averaging buffer and provide an immediate update on the display as necessary. This is especially useful when measuring samples with large differences in turbidity.

Maintenance

⚠ DANGER

Multiple hazards. Only qualified personnel must conduct the tasks described in this section of the document.

Clean the instrument

Keep the instrument clean to get continuous and accurate operation.

NOTICE

Never use cleaning agents such as turpentine, acetone or similar products to clean the instrument including the keypad.

1. Turn the instrument off and disconnect the power cord.
2. Clean the surface of the instrument with a soft, moist cloth and a weak soap solution.
3. Dry the surface of the instrument with a lint-free cloth.

Replace the LED light source

The light source, light emitting diode (LED), is not user replaceable. Contact Customer Service for LED replacement.

Replace a fuse

⚠ DANGER


Fire hazard. Use the same type and current rating to replace fuses.

Replacement parts:

- Fuse for 115 V operation, time-delay, 250 V, 1.6 A (3030700), or
- Fuse for 230 V operation, time-delay, 250 V, 1.6 A (3030600)

To replace a fuse, refer to the illustrated steps in [Figure 6](#).

Figure 6 Replace a fuse

Troubleshooting

Refer to the tables in this section for error codes, diagnostic codes, common problem messages or symptoms, possible causes and corrective actions.

Error codes

[Table 3](#) lists the error codes shown for different conditions. Error codes identify instrument malfunction or operator error.

The instrument continues operation in an error condition.

Push **ENTER** to clear an error code from the display.

Note: Any calibration being calculated when an error occurs, is discarded. The old calibration is kept.

Table 3 Error codes

Error	Description	Solution
ERR01	The turbidity of the dilution water is greater than 0.5 NTU.	Start the calibration again with lower turbidity dilution water. Note: Ignore ERR01 when the sample cell diameter is less than 25 mm. Push UNITS/Exit to go back to measurement mode.
ERR02	- Two calibration standards have the same value. - The difference between two calibration standards is less than 60.0 NTU. - The turbidity of Standard 1 is too low (less than 10 NTU).	- Inspect the preparation of standards. - Do the calibration again. Note: Ignore ERR02 when the sample cell diameter is less than 25 mm. Push UNITS/Exit to go back to measurement mode.
ERR03	Low light error	- Put the sample in the instrument again. - Make sure that the lamp light is on. - Make sure that an object is not in the light path. - Do sample dilution if necessary.
ERR04	Memory malfunction	- Turn the instrument off and then back on. - Contact Technical Support if the error occurs again.
ERR05	A/D is over the range	- Make sure that the light shield is closed. - Contact Customer Service if necessary.

Table 3 Error codes (continued)

Error	Description	Solution
ERR06	A/D is under the range	- Make sure that no object is in the light path. - Contact Customer Service if necessary.
ERR07	Light leak	- Make sure that the cover for the sample cell compartment is closed. - Turn the instrument off and then back on.
ERR09	Printer time out error	- Make sure that the external printer is connected correctly. - Make sure that the external printer is selected (online).
ERR10	System voltage out of range	- Turn the instrument off and then back on. - Contact Customer Service if the error occurs again.
ERR11	System loop test error	- Turn the instrument off and then back on. - Contact Customer Service if the error occurs again.

Diagnostic codes

Table 4 lists the diagnostic codes that are used to get information about instrument operation when instrument operation is in doubt.

To do a diagnostic test:

1. Push **SETUP**.
2. Use the arrow keys to enter a diagnostic code.
3. Push **ENTER** to show the diagnostic value.
4. Push **UNITS/Exit** to go back to Measurement mode.

Note: To print a diagnostic report, hold down **PRINT**, then turn the instrument on.

Table 4 Diagnostic codes

Code	Display	Description
21	Pr In	Printer test
22	Test results are shown.	Display test
23	Test results are shown.	Keyboard test
24	Test results are shown.	Memory test

Delete calibration data

To delete any calibration data entered by the user:

1. Turn off the instrument.
2. Push and hold **CAL**.
3. Turn on the instrument.

The CAL? light flashes. The instrument starts in Calibration mode.

4. Calibrate the instrument before use.

Flashing 9s

When manual ranging is selected, the display will flash all 9s when the sample being measured is greater than the selected range.

When automatic ranging is selected, the display will flash 9s when the sample is greater than the maximum range of the instrument.

Flashing 0s

When manual ranging is selected, the display will flash all 0s when the sample measured is less than the selected range.

When automatic ranging is selected, the display will flash all 0s when the measurement is less than the range of the instrument or a negative value. Calibrate the instrument.

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P.O. Box 389, Loveland, CO 80539-0389 U.S.A.
(800) 227-4224 (U.S.A. only)