

modified colitag™

Presence-Absence



Get EPA-approved
Presence/Absence
results faster.

DW

Product Overview

Hach Modified Colitag™ is the fastest presence-absence test in the market—with a 16- or 22-hour testing timeframe. This allows drinking water facilities to begin a test at 5:00 p.m. and get EPA-approved results by 9:00 am the following morning.

The product also provides unparalleled flexibility by offering two testing parameters in one product: a 16-hour or a 22-hour test. Being able to select between two tests in one product, allows drinking water facilities to set up the best timeframes for their staffing and to adjust tests as needed, without changing procedures. In addition, Hach's Modified Colitag provides users with valid results for up to 32 hours after the completion of the test. This best-in-class testing timeframe allows for both flexible workflow and testing over a weekend.

E

Features & Benefits

The fastest on the market at 16 hours

- Fastest EPA-approved test for coliforms and E. coli
- Test as late as 5:00 PM and still have your results by 9:00 AM the next day

One product that supports 16- or 22-hour test methods

- Allows for a "standard" test that requires less set-up, or a "rapid" test when necessary
- One media for both tests reduces complexity in certification and inventory

Results are valid for 32 hours

- Saturday's tests can be read on Monday—no need to work all weekend
- Results stay valid for 800% longer than the competition

FB

DW = drinking water WW = wastewater municipal PW = pure water / power
IW = industrial water E = environmental C = collections FB = food and beverage



Be Right™

Hach Modified Colitag Procedure

Step 1: Snap to open Hach Modified Colitag Pack



Step 2: Add Hach Modified Colitag to 100 mL water sample



Step 3: Shake sample, media will not dissolve immediately



Step 4: Incubate for 16-22 Hours at 35°C ±0.5°C



Step 5: Check for yellow color, indicating coliforms



Step 6: Check for fluorescence with UV light, indicating *E. coli*



Ordering Information

- 8888001** Comparitor for use with Hach Modified Colitag
8888004 100 Pack, Snap-Pack
8888005 200 Pack, Snap-Pack
8888006 120 mL Sterile bottles, 100 Pack
2184300 UV Lamp, 365 nm

Lit. No. 2715 Rev 1

A12 Printed in U.S.A.

©Hach Company, 2012 All rights reserved.

In the interest of improving and updating its equipment, Hach Company reserves the right to alter specifications to equipment at any time.

At Hach, it's about learning from our customers and providing the right answers. It's more than ensuring the quality of water—it's about ensuring the quality of life. When it comes to the things that touch our lives...

Keep it pure.

Make it simple.

Be right.

For current price information, technical support, and ordering assistance, contact the Hach office or distributor serving your area.

In the United States, contact:

HACH COMPANY World Headquarters
P.O. Box 389
Loveland, Colorado 80539-0389
U.S.A.
Telephone: 800-227-4224
Fax: 970-669-2932
E-mail: orders@hach.com
www.hach.com

U.S. exporters and customers in Canada, Latin America, sub-Saharan Africa, Asia, and Australia/New Zealand, contact:

HACH COMPANY World Headquarters
P.O. Box 389
Loveland, Colorado 80539-0389
U.S.A.
Telephone: 970-669-3050
Fax: 970-461-3939
E-mail: intl@hach.com
www.hach.com

In Europe, the Middle East, and Mediterranean Africa, contact:

HACH LANGE GmbH
Willstätterstraße 11
D-40549 Düsseldorf
GERMANY
Tel: +49 (0) 211 5288-0
Fax: +49 (0) 211 5288-143
E-mail: info@hach-lange.de
www.hach-lange.com



Be Right™